

Sensing and Sensomics

*the chemistry behind olfactory
pleasure*



Chiara Cordero

Food Chemistry and Pharmaceutical Biology
Dipartimento di Scienza e Tecnologia del Farmaco
Università degli Studi di Torino

Overview



Foreword

Sensory Quality - Food Selection - Quality of life

The Olfactory Sense

- ✓ Olfactory Receptors and Receptor Code
- ✓ The human Odorome

Modern “*omics*” in food research: *sensomics*

Analytical strategies designed to exploit the complexity

Advanced *sensomics* platforms

- ✓ Multidimensionality as a key-factor
- ✓ Unrevealing the Chemical Odor Code of Milk




Concluding remarks

Behind pleasure, a further key for well-being

Foreword



Plant volatilome



Food odorome




"...it has become possible to include the chemical space in system biology, i.e. systems chemical biology."

Oprea TI, Tropsha A, Faulon JL, Rintoul MD (2007) Systems chemical biology. Nat Chem Biol 3: 447–450

Sensory Quality and Food Selection and much more




Daily life problems in patients with smell disorders





Complex and intriguing task:

1. Native flavor compounds trace level
2. Olfactory perception is triggered by ligand-receptor interactions
3. One ligand interacts with more than one receptor
4. Receptors activation generates a complex patter of signals




The Olfactory Sense



olfactory epithelium

STEP 3



olfactory bulb

STEP 4

Chemistry behind olfaction:

- ✓ Identification of ligands for olfactory receptors (pharmacological space)
- ✓ Identification of chemical features related odor quality and odor threshold (SAR studies – medicinal chemistry space)
- ✓ Characterization of perireceptor events that modulate odorant-receptor interaction (enzymatic activity)
- ✓ Identification of chemical patterns related with the combinatorial code of olfaction



Olfactory receptors



The Nobel Prize Medal for Physiology or Medicine Year 2004



Richard Axel




Linda B. Buck


"for their discoveries of odorant receptors and the organization of the olfactory system"

ORs are members of a superfamily of up to 1000 different **G-protein coupled receptors (G_{olf})** with seven trans-membrane (7TM) domains.

The location of odorant binding is thought to be a hydrophobic pocket in trans-membrane regions 3, 4, and 5 of the seven-membrane-spanning receptor.



Canonical pathway of signal transduction in olfactory sensory neurons (OSN)



Ligand-Receptor interaction

Dose response curves

Functional expression of ORs in heterologous expression systems, i.e., human embryonic kidney cells (HEK293).

Acetophenone (Acp) and **Benzaldehyde** (Bnz)
Agonists (different sensitivity -10 fold dose shift)
at M71 OR

M71-expressing olfactory sensory neurons (OSNs) are evidenced by green fluorescence.




Figure 3 Recording odorant-induced responses from olfactory sensory neurons expressing the M71 OR. **(A)** The introduction of IRES-tauGFP into the mouse genome following the coding sequence for the M71 OR allowed the identification and visualization of all M71-expressing olfactory sensory neurons, which are found in the most dorsal zone on the olfactory turbinates. Scalebar, 250 μm . **(B)** Screening of GFP-positive olfactory sensory neurons using Ca^{2+} imaging and a variety of odorants revealed the M71 agonists acetophenone (Acp) and benzaldehyde (Bnz). Note the ~ 10 -fold shift difference in sensitivity. Modified from Bozza et al. (2002) with permission from the Society of Neuroscience.



Ligand-Receptor interaction "Odotope theory"

The Olfactory Sense
Odotope-theory and
concentration thresholds




A combination of activated receptors is responsible for any one smell – **Receptor Code**


Receptors recognize small structural features on each molecule and the brain is responsible for processing the combined signal into an interpreted smell



Caramel-like




Furaneol




Maltol

Spicy flavor



Sotolon



Ethyl-Sotolon

Ligand-Receptor interaction “Odor potency”

The Olfactory Sense
Odotope-theory and
concentration thresholds



JOURNAL OF
AGRICULTURAL AND
FOOD CHEMISTRY

Article
pubs.acs.org/JAFC

Influence of the Chemical Structure on Odor Qualities and Odor Thresholds in Homologous Series of Alka-1,5-dien-3-ones, Alk-1-en-3-ones, Alka-1,5-dien-3-ols, and Alk-1-en-3-ols

Katja Lorber,^{#,†} Peter Schieberle,[§] and Andrea Buettner^{*,#,\dagger}


odorant	odor quality ^b	threshold in air ^a (ng/L)	
		by GC-O	data from lit. ^c
pent-1-en-3-one	pungent, fish-like 	0.17	20 ²²
hex-1-en-3-one	metallic	0.02	0.02–0.04, ²³ 0.7 ²²
hept-1-en-3-one	vegetable-like, metallic	0.08	0.7 ²²
oct-1-en-3-one	mushroom-like, metallic 	0.04	0.03–1.12, ^{24,25} 0.1 ²²
non-1-en-3-one	mushroom-like 	0.008	0.008, ⁴⁰ 0.02 ²²
dec-1-en-2-one	mushroom-, geranium-, herb-like 	0.08	1 ²²
undec-1-en-3-one	herb-, citrus-like, soapy 	0.15	1 ²²
dodec-1-en-3-one	herb-, citrus-like, soapy 	33	180 ²²
tridec-1-en-3-one	citrus-like, soapy 	166	nr
tetradec-1-en-3-one	soapy, plastic-like 	771	nr

Higher level of information is needed to correctly interpret the complex phenomena behind olfactory perception.

1. Humans ≈ 400 different ORs



2. Each OSN expresses only one kind ORs




Selectivity study on 60 olfactory neurons tested with 20 different odorants. Each neuron is equivalent to one column, each odorant to one row. The size of the dots symbolizes the number of action potentials.

Neurons 26 and 58 react specifically on cineole and camphor.


Neuron 4 reacts on all odorants except isoamyl acetate.



3. ORs are not highly specific multiple activation/modulation



Receptor Code
Complex activity pattern of many differently tuned olfactory neurons that project in specific regions of the olfactory bulb



Olfactory Bulb Odor Map

Ligand: **Isoamyl acetate**

Species: Rat – Strain: Sprague-Dawley


Exposure time: 128 sec

Brain Region: Olfactory bulb



Sub Region: Glomerular layer

Side: Right

View: Lateral



Topography of
the activation
pattern



To lateral
olfactory tract
CNS integration



Functional MRI images revealing the olfactory bulb activation pattern for ISOAc at different concentrations (4500 – 520 – 41 ppm)

Modeling of the odorant human combinatorial coding

OPEN ACCESS Freely available online

PLOS ONE



The Olfactory Sense
Bioinformatics



Identification of Odorant-Receptor Interactions by Global Mapping of the Human Odorome

Karine Audouze¹, Anne Tromelin², Anne Marie Le Bon², Christine Belloir², Rasmus Koefoed Petersen³, Karsten Kristiansen³, Søren Brunak¹, Olivier Taboureau^{1,4*}

April 2014



Node size: number of odorants known to bind to the OR (83 ORs and 938 associations).

Width of the edges: weighted score is the strength of the link between two ORs based on the number of shared compounds.

Mapping of odor descriptions on the human odorome by an association score based on the number of compounds associated to an odor.



Q. Starting from food composition, how can we tackle (at a molecular basis) these complex phenomena, how can we objectify odor perception?
A. Comprehensive chemical investigation at multiple level of information

Foodomics¹

investigates food and nutrition domains (food science and nutrition) combining information from emerging working areas, **advanced analytical techniques** and bioinformatics

Sensomics² - Flavoromics^{3,4} - Flavour metabolomics⁵

correlate **analytical data to perception** extend investigation to all possible *stimuli* of a multimodal perception **collecting data on volatiles, semi-volatiles** and non-volatiles **in an efficient manner**

[1] Herrero M. et al (2012) Mass Spectrom Rev. 31(1):49-69

[2] Schieberle, P., Hofmann, T. (2003) Chemie in Unserer Zeit, 37 (6), pp. 388-401.

[3] G. A. Reineccius. (2008) In Abstracts of Papers, 235th American Chemical Society, New Orleans, LA, USA, pp. AGFD-061.

[4] Charve J et al (2011) Flavour Fragrance J. 26(6):429-440.

[5] de Vos C. H. R. et al (2008) Proceedings of the 12th Weurman Symposium, Interlaken, Switzerland pp. 573-580

Collaboration with Prof. Peter Schieberle
Food Chemistry Department
Technical University of Munich





Targeted analysis

Chemical features are just known-markers and analytical efforts are directed to few informative analytes

informative potential limited to known attributes

Information from unknowns is lost and untreated for routine practice

Classic Sensomics
Focus on few relevant odor active compounds



Untargeted analysis

Each and every sample component is a potential *informative chemical feature* all detectable chemical information is considered

informative potential is extended to the whole chemical pattern

Exploits the hidden information enveloped in the sample fingerprint

Deep insights into the complex phenomena of sensory perception
Detailed chemical characterization



Advanced *Sensomics* approaches aims at merging discrete and time-consuming operations in one step by adopting multidimensional analytical platforms

1. extraction/isolation of volatile fraction compounds
2. concentration of extracts
3. pre-separation of volatile compounds
4. main separation of volatile compounds
5. localization of intense aroma compounds
6. identification of aroma compounds
7. quantification of aroma compounds
8. evaluation of their aroma contribution
9. aroma recombination

Sensu lato multidimensional techniques:

Liquid Chromatography coupled with high/low resolution Mass Spectrometry (LC-MSⁿ)

Gas Chromatography coupled with high/low resolution Mass Spectrometry (GC-MS)

Sensu stricto multidimensional techniques:

Two-Dimensional Comprehensive GC (GC × GC) coupled with MS and Tandem MS

Two-Dimensional Comprehensive LC (LC × LC) coupled with MS and Tandem MS

Omics Analytical Platforms

LC-MS LC-MSⁿ LC×LC-MS LC×LC-MSⁿ
HS-MS GC-MS GC×GC-MS GC×GC-MSⁿ

Chemical information

Identity

exact mass assignment
fragmentation pattern
multiple reaction monitoring


Physicochemical behavior

volatility, polarity, partition coefficient,
solubility

Mass quantitation

absolute concentration and relative
abundance

Bio-guided screening Olfactometry & Sensory analysis



Aroma, taste, texture, aspect, etc...
orthogonal and **complementary**
treated together with chemical data





interpret chemical speciation
Identifying the Chemical Odor Code






Omics Analytical Platforms
 LC-MS LC-MSⁿ LC×LC-MS LC×LC-MSⁿ
 HS-MS GC-MS **GC×GC-MS** **GC×GC-MSⁿ**


1st dimension: **Two-dimensional comprehensive gas chromatography**



Second dimension: polyacrylamide gel electrophoresis (PAGE)




Further information dimensions
 Molecular identity (MS spectrum)
 Odor descriptor (Olfactometry)
Multidimensional information space



Omic Analytical Platforms

LC-MS LC-MSⁿ LC×LC-MS LC×LC-MSⁿ
HS-MS GC-MS **GC×GC-MS** **GC×GC-MSⁿ**

Modern omics disciplines
Advanced Sensomics




Aldehydes

2-Methyl Butanal
Crotonaldehyde
Capronaldehyde
2E-Heptenal
Benzaldehyde
3-Carboxaldehyde Thiophene
Caprylaldehyde
Phenyl Acetaldehyde
Nonanal
2-Nonenal

Omic Analytical Platforms

LC-MS LC-MSⁿ LC×LC-MS LC×LC-MSⁿ
HS-MS GC-MS **GC×GC-MS** **GC×GC-MSⁿ**

Modern omics disciplines
Advanced Sensomics




Furans

2-Methylfuran
THF
Dihydro-2(3H)-Furanone
2-Furancarboxaldehyde
2-Furanemethanol
2-Acetylfuran
2(5H)-FURANONE
5-Methyl-2-Furanemethanol
5-Methyl-2-Furanecarboxaldehyde
2,5-Dimethyl-3(2H)Furanone

Omic Analytical Platforms

LC-MS LC-MSⁿ LC×LC-MS LC×LC-MSⁿ
HS-MS GC-MS GC×GC-MS GC×GC-MSⁿ

Modern omics disciplines
Advanced Sensomics




Acids

Formic acid
Acetic acid
Propionic acid
Butyric acid
Valeric acid
Isovaleric acid
2-methyl butanoic acid
Capronic acid
Enanthic acid
Caprylic acid

Omic Analytical Platforms

LC-MS LC-MSⁿ LC×LC-MS LC×LC-MSⁿ
 HS-MS GC-MS **GC×GC-MS** **GC×GC-MSⁿ**

Modern omics disciplines
 Advanced Sensomics



Pyrazines

Pyrazine
 Methylpyrazine
 2,5-Dimethylpyrazine
 2,6-Dimethylpyrazine
 2-Ethylpyrazine
 2,3-Dimethylpyrazine
 2-Ethyl-6-Methylpyrazine
 2,3,5-Trimethylpyrazine
 2-Ethyl-5-Methylpyrazine
 2-Ethyl-3-Methylpyrazine
 2-Ethenyl-5-methylpyrazine
 2-Ethenyl-6-methylpyrazine
 2-Acetyl pyrazine
 2-Ethyl-3,6-dimethylpyrazine
 2,3-Diethylpyrazine
 3-Ethyl-2,5-diMethyl pyrazine
 2-Isopropenylpyrazine
 2-Acetyl-6-methylpyrazine
 2-Acetyl-5-methylpyrazine
 2,3-Diethyl-5-methylpyrazine
 (E)-2-Methyl-5-(1-propenyl)pyr
 2-Acetyl-3-ethylpyrazine

Biometric Fingerprinting¹


Modern omics disciplines
Advanced Sensomics



Each constituent (chemical feature) has a distinct position over the 2D plane and a unique information pattern (MS spectrum) that enables **identification and quantitation**

Comparison with stored templates

Two-dimensional separation pattern (2D chromatogram) of hazelnut samples obtained by two-dimensional comprehensive gas chromatography- mass spectrometry (GCxGC-MS)



This 2D pattern is unique and informative, collects data on multiple sample properties, including **odorants distribution.**

The pattern is the sample

(Detailed) Profiling²

Identification

Rank	NC0041942	Lat	LCF
2)	NC0965119A	2335	F#6
3)	NC0843399A	2315	F#6
4)	NC0889984A	2190	F#6
5)	NC0172646A	2170	F#6
6)	NC0948401A	2160	F#6
7)	NC0163889A	2080	F#6
8)	NC0199511A	2075	F#6
9)	NC0814921A	2075	F#6
10)	NC0062100G	2045	F#6
11)	NC0339094A	2040	F#6
12)	NC0529094A	2025	F#6
13)	NC0782084A	2020	F#6
14)	NC0058654H	2010	F#6

Case Information
Case ID : N000001221
Resp ID : NC0771078A
Resp : 18
Lists : 1

Disposition Case Cancel Verification

Case Selection Control
Status : Case 1 of 1



Filter: N00000122
Query: local
Expand

Sensomics advanced Platform
Sample preparation Separation-Detection
Sensory screening Data elaboration

Modern omics disciplines
 Advanced Sensomics



Collaboration with
Prof. S.E. Reichenbach
Computer Science and Engineering
Department University of Nebraska



Data Acquisition: Agilent Chem Station ver D.02
Data elaboration : GC-Image ver 2.1b4

Agilent 6890-5975C MS
 operating in EI mode (70 eV) Scan range 35/240 m/z, Scan rate 10,000 amu/s

Zoex KT2004 loop modulator
 Cryogenic: liquid nitrogen, Massflow ctr. Optimode™ by SRA Italy

GC-O sniffing port:
 design by Grosh et al at FID outlet before ²D column.

GC×GC-MS: injector temperature: 250°C, injection mode: split, ratio: 1/10; carrier gas: helium, initial head pressure 260 KPa;
Column set: ¹D Sol-GelWax® (30 m, 0.25 mm ID, 0.25 µm d_f) (SGE Analytical Science, Ringwood, Victoria, Australia) and ²D OV1701 (1 m, 0.1 mm ID, 0.1 µm d_f) (Mega, Legnano (Milan), Italy); loop capillary 1 m, 0.1 mm ID, 0.1 µm d_f
GC oven: 40°C (1 min) to 170°C at 10°C/min and to 260°C at 50°C/min (1 min)

Cordero, C., et al. "Targeted and Non-Targeted Approaches for Complex Natural Sample Profiling by GC×GC-qMS." *Journal of chromatographic science* 48.4 (2010): 251-61.

Cordero, C., et al. "Profiling Food Volatiles by Comprehensive Two-Dimensional Gas chromatography Coupled with Mass Spectrometry: Advanced Fingerprinting Approaches for Comparative Analysis of the Volatile Fraction of Roasted Hazelnuts (*Corylus Avellana* L.) from Different Origins." *Journal of Chromatography A* 1217.37 (2010): 5848-58.

Kiefl, J., et al. "Performance Evaluation of Non-Targeted Peak-Based Cross-Sample Analysis for Comprehensive Two-Dimensional Gas Chromatography-Mass Spectrometry Data and Application to Processed Hazelnut Profiling." *Journal of Chromatography A* 1243 (2012): 81-90.

Reichenbach, S. E., et al. "Features for Non-Targeted Cross-Sample Analysis with Comprehensive Two-Dimensional Chromatography." *Journal of Chromatography A* 1226 (2012): 140-8.

Reichenbach, S. E., et al. "Reliable Peak Selection for Multisample Analysis with Comprehensive Two-Dimensional Chromatography." *Analytical Chemistry* 85.10 (2013): 4974-81.

Whole and nonfat dry milk: a case study



Maillard reaction products (including sulphur derivatives)

Strecker aldehydes

Lipid oxidation products (enzymatic and ozonolysis during spray-drying)

Products of TG hydrolysis (free fatty acids, including short chain fatty acids)

Lactones and furanones

Aromatic AA degradation (phenolic derivatives responsible of cowy and fecal flavors)

Terpenoids and norisoprenoids


Target chemical features

Volatility range (25°C) from 49.3 of 3-methylbutanal to 1.72E-5 Torr of (*E*)-2-octadecenal

Log P values range: 0.32 of acetic acid to 8.05 of octadecanal

Water solubility (pH 6 – 25°C) range: from 16.6 of acetic acid to 5.6E-06 mol/L of 2-octadecanone

Odor thresholds (orthonasal in water at 25°C) 0.15 ng/g for the (*E*)-2-nonenal



Sensory features

odor descriptors

off-flavors (shelf life)

positive attributes responsible of the aroma



Highly volatiles - Orthonasal perception

Static Headspace
 Temperature - Time - Phase ratio (β)

Dynamic Headspace
 breakthrough volume

S-HS	HS-SPME Selective extraction	HSSE Larger volumes	HS-MMSE High surface
------	---------------------------------	------------------------	-------------------------

D-HS
 Larger sampling capacity

Sorption PDMS PA	Adsorption CW	Mixed mode	Sorption PDMS	Mixed mode PDMS/ Carbopack	Adsorption ODS Carbopack
------------------------	------------------	---------------	------------------	-------------------------------------	--------------------------------

Sorption PDMS
 (foams and particles)

Adsorption
 (Tenax TA)

SPME	SBSE	MMSE
------	------	------

In-solution sampling
 Temperature - Time - Phase ratio (β)

Evaporation in oral cavity - Retronasal perception




S-HS




Sampling: 50 (°C)
Syringe temperature: 60 (°C) Incubation: 40 min
Fill speed: 40 (mL/min)
Sample volume: 1 (mL) Fill volume: 2 (ml)
Injection speed: 100 (mL/min)
Flush time: 2 (min)
System HT2100H (HTA srl - Brescia Italy)

HS-SPME
Selective extraction



Sampling temperature: 50°C
Sampling time: 40 min
Vial volume: 20 mL
Fiber PDMS 100 µm
Polyacrylate 85 µm
Carbowax 60 µm
DVB/CAR/PDMS 50/30 µm 2 cm
Supelco Bellefonte


HSSE
Larger volumes



Sampling temperature: 50°C
Sampling time: 40 min
Vial volume: 20 mL


Conventional PDMS Twister (length: 20 mm, thickness: 0.5 mm).
Dual phase Twisters PDMS-Carbopack B (length: 20 mm, thickness: 0.5 mm)
Gerstel Mulheim DE.

HS-MMSE
High surface



Sampling temperature: 50°C
Sampling time: 40 min
Vial volume: 20 mL
Monotrap Rod type: 2.9 mmOD x 5 mm height
Octadecylsilica-Active Carbon Black GL
Sciences Inc. Tokyo Japan



D-HS
Larger sampling capacity



Sampling temperature: 40°C
Trap temperature: 25°C Time: 20 min
Stripping gas: N₂ at 10 mL/min
Vial volume: 20 mL
PDMS particles and foam by Gerstel Mulheim DE
Tenax TA (60/80 mesh) Supelco Bellefonte

Static Headspace (8)

Dynamic Headspace (2)



SPME
Sorption

SPME
Adsorption

SPME
Mixed
mode


HSSE
Sorption

HSSE
Mixed
mode

MMSE
Adsorption

Sorption PDMS

Adsorption



In-solution sampling (6)

Whole and nonfat milk (2)
Sampling replicates (3)
 $8 \times 2 \times 6 \times 3 = 96$ 2D patterns




Bio-guided Chromatographic fingerprinting
Detailed untargeted Profiling
Chemical Odor Code of Milk





Bio-guided fingerprinting **Bio-guided screening of samples pattern**

Low resolution fingerprinting
based on *region features*¹




specific regions of the chromatographic plane where detector
(single or multi channel) intensities are summed together
Information is extracted , aligned and compared across samples



these regions can be fixed *a priori* (grids)
enables bio-guided screening of the most informative regions of
the chromatographic retention space


[1] Reichenbach S.E., Tian X., Cordero C., Tao T. (2012) J. Chromat. A, 1226 :140– 148

Region features




Bio-guided Chromatographic Fingerprinting¹


1. comprehensive subdivision of the chromatographic plane



- ✓ uniform grid with a fixed number of contiguous and non-overlapping panels corresponding to 1D LRI
- ✓ GC(O)×GC-MS screening



Picture from:
TrAC (2011) 30(11):1756–1770



Detailed untargeted Profiling Comprehensive Template Matching Fingerprinting¹


Peak feature

To improve the specificity of the process one has to limit correspondences to peaks resulting from the same analyte




A	B	C	D	E	F	
BlobID	Compound ID	Peak I (min)	Peak II (sec)	Peak Volume	LRI	Reference MS
374	121	19.02	3.79	245302	975	35.0,2.0;36.0,80.0;37.0,115.0;38.0,38.0;39.0,12208.0;40.0,837.0;41.0,25786.0;42.0,7820.0;43.0,1505
332	102	19.02	1.52	110655	975	36.0,40.0;37.0,157.0;38.0,183.0;39.0,4372.0;40.0,894.0;41.0,4512.0;42.0,5043.0;43.0,12806.0;44.0,6
172	108	19.09	2.74	425781	976	35.0,13.0;36.0,24.0;37.0,1695.0;38.0,1986.0;39.0,9072.0;40.0,291.0;41.0,5354.0;43.0,152234.0;44.0,
245	106	19.15	2.23	341990	977	36.0,15.0;37.0,507.0;39.0,42339.0;40.0,2480.0;41.0,11773.0;42.0,1162.0;43.0,17763.0;44.0,651.0;45
77	107	19.29	2.65	2796671	980	35.0,680.0;36.0,43.0;37.0,9139.0;38.0,20531.0;39.0,110483.0;40.0,15765.0;41.0,83686.0;42.0,14859
90	103	19.29	0.93	4295315	980	35.0,2.0;36.0,203.0;37.0,58686.0;38.0,122582.0;39.0,380440.0;40.0,173634.0;41.0,21258.0;42.0,188
456	120	19.35	3.66	92337	981	36.0,45.0;37.0,561.0;38.0,393.0;39.0,12797.0;40.0,1557.0;41.0,25775.0;42.0,1061.0;43.0,1658.0;44.0,
292	105	19.42	1.81	152393	982	36.0,4.0;37.0,19.0;38.0,123.0;39.0,5130.0;40.0,478.0;41.0,5227.0;42.0,5363.0;43.0,6801.0;44.0,380.0
343	110	19.55	2.69	218605	985	35.0,81.0;36.0,79.0;37.0,451.0;38.0,1245.0;39.0,16581.0;40.0,501.0;41.0,5990.0;42.0,118.0;43.0,340
252	119	19.82	2.82	320953	989	36.0,161.0;37.0,750.0;38.0,1302.0;39.0,27531.0;40.0,4668.0;41.0,18879.0;42.0,6904.0;43.0,39243.0;
1	112	19.89	2.53	7.31E+07	990	35.0,34.0;36.0,250.0;37.0,111415.0;38.0,373491.0;39.0,1511619.0;40.0,1917573.0;41.0,297586.0;42.
127	104	19.89	1.05	696840	990	35.0,24.0;36.0,109.0;37.0,1720.0;38.0,601.0;39.0,12091.0;40.0,2215.0;41.0,1907.0;42.0,3782.0;43.0,
247	117	20.09	2.74	703966	994	36.0,37.0;37.0,12.0;38.0,12.0;39.0,12.0;40.0,12.0;41.0,12.0;42.0,12.0;43.0,12.0;44.0,12.0;45.0,12.0;46.0,12.0;47.0,12.0;48.0,12.0;49.0,12.0;50.0,12.0;51.0,12.0;52.0,12.0;53.0,12.0;54.0,12.0;55.0,12.0;56.0,12.0;57.0,12.0;58.0,12.0;59.0,12.0;60.0,12.0;61.0,12.0;62.0,12.0;63.0,12.0;64.0,12.0;65.0,12.0;66.0,12.0;67.0,12.0;68.0,12.0;69.0,12.0;70.0,12.0;71.0,12.0;72.0,12.0;73.0,12.0;74.0,12.0;75.0,12.0;76.0,12.0;77.0,12.0;78.0,12.0;79.0,12.0;80.0,12.0;81.0,12.0;82.0,12.0;83.0,12.0;84.0,12.0;85.0,12.0;86.0,12.0;87.0,12.0;88.0,12.0;89.0,12.0;90.0,12.0;91.0,12.0;92.0,12.0;93.0,12.0;94.0,12.0;95.0,12.0;96.0,12.0;97.0,12.0;98.0,12.0;99.0,12.0;100.0,12.0;101.0,12.0;102.0,12.0;103.0,12.0;104.0,12.0;105.0,12.0;106.0,12.0;107.0,12.0;108.0,12.0;109.0,12.0;110.0,12.0;111.0,12.0;112.0,12.0;113.0,12.0;114.0,12.0;115.0,12.0;116.0,12.0;117.0,12.0;118.0,12.0;119.0,12.0;120.0,12.0;121.0,12.0;122.0,12.0;123.0,12.0;124.0,12.0;125.0,12.0;126.0,12.0;127.0,12.0;128.0,12.0;129.0,12.0;130.0,12.0;131.0,12.0;132.0,12.0;133.0,12.0;134.0,12.0;135.0,12.0;136.0,12.0;137.0,12.0;138.0,12.0;139.0,12.0;140.0,12.0;141.0,12.0;142.0,12.0;143.0,12.0;144.0,12.0;145.0,12.0;146.0,12.0;147.0,12.0;148.0,12.0;149.0,12.0;150.0,12.0;151.0,12.0;152.0,12.0;153.0,12.0;154.0,12.0;155.0,12.0;156.0,12.0;157.0,12.0;158.0,12.0;159.0,12.0;160.0,12.0;161.0,12.0;162.0,12.0;163.0,12.0;164.0,12.0;165.0,12.0;166.0,12.0;167.0,12.0;168.0,12.0;169.0,12.0;170.0,12.0;171.0,12.0;172.0,12.0;173.0,12.0;174.0,12.0;175.0,12.0;176.0,12.0;177.0,12.0;178.0,12.0;179.0,12.0;180.0,12.0;181.0,12.0;182.0,12.0;183.0,12.0;184.0,12.0;185.0,12.0;186.0,12.0;187.0,12.0;188.0,12.0;189.0,12.0;190.0,12.0;191.0,12.0;192.0,12.0;193.0,12.0;194.0,12.0;195.0,12.0;196.0,12.0;197.0,12.0;198.0,12.0;199.0,12.0;200.0,12.0;201.0,12.0;202.0,12.0;203.0,12.0;204.0,12.0;205.0,12.0;206.0,12.0;207.0,12.0;208.0,12.0;209.0,12.0;210.0,12.0;211.0,12.0;212.0,12.0;213.0,12.0;214.0,12.0;215.0,12.0;216.0,12.0;217.0,12.0;218.0,12.0;219.0,12.0;220.0,12.0;221.0,12.0;222.0,12.0;223.0,12.0;224.0,12.0;225.0,12.0;226.0,12.0;227.0,12.0;228.0,12.0;229.0,12.0;230.0,12.0;231.0,12.0;232.0,12.0;233.0,12.0;234.0,12.0;235.0,12.0;236.0,12.0;237.0,12.0;238.0,12.0;239.0,12.0;240.0,12.0;241.0,12.0;242.0,12.0;243.0,12.0;244.0,12.0;245.0,12.0;246.0,12.0;247.0,12.0;248.0,12.0;249.0,12.0;250.0,12.0;251.0,12.0;252.0,12.0;253.0,12.0;254.0,12.0;255.0,12.0;256.0,12.0;257.0,12.0;258.0,12.0;259.0,12.0;260.0,12.0;261.0,12.0;262.0,12.0;263.0,12.0;264.0,12.0;265.0,12.0;266.0,12.0;267.0,12.0;268.0,12.0;269.0,12.0;270.0,12.0;271.0,12.0;272.0,12.0;273.0,12.0;274.0,12.0;275.0,12.0;276.0,12.0;277.0,12.0;278.0,12.0;279.0,12.0;280.0,12.0;281.0,12.0;282.0,12.0;283.0,12.0;284.0,12.0;285.0,12.0;286.0,12.0;287.0,12.0;288.0,12.0;289.0,12.0;290.0,12.0;291.0,12.0;292.0,12.0;293.0,12.0;294.0,12.0;295.0,12.0;296.0,12.0;297.0,12.0;298.0,12.0;299.0,12.0;300.0,12.0;301.0,12.0;302.0,12.0;303.0,12.0;304.0,12.0;305.0,12.0;306.0,12.0;307.0,12.0;308.0,12.0;309.0,12.0;310.0,12.0;311.0,12.0;312.0,12.0;313.0,12.0;314.0,12.0;315.0,12.0;316.0,12.0;317.0,12.0;318.0,12.0;319.0,12.0;320.0,12.0;321.0,12.0;322.0,12.0;323.0,12.0;324.0,12.0;325.0,12.0;326.0,12.0;327.0,12.0;328.0,12.0;329.0,12.0;330.0,12.0;331.0,12.0;332.0,12.0;333.0,12.0;334.0,12.0;335.0,12.0;336.0,12.0;337.0,12.0;338.0,12.0;339.0,12.0;340.0,12.0;341.0,12.0;342.0,12.0;343.0,12.0;344.0,12.0;345.0,12.0;346.0,12.0;347.0,12.0;348.0,12.0;349.0,12.0;350.0,12.0;351.0,12.0;352.0,12.0;353.0,12.0;354.0,12.0;355.0,12.0;356.0,12.0;357.0,12.0;358.0,12.0;359.0,12.0;360.0,12.0;361.0,12.0;362.0,12.0;363.0,12.0;364.0,12.0;365.0,12.0;366.0,12.0;367.0,12.0;368.0,12.0;369.0,12.0;370.0,12.0;371.0,12.0;372.0,12.0;373.0,12.0;374.0,12.0;375.0,12.0;376.0,12.0;377.0,12.0;378.0,12.0;379.0,12.0;380.0,12.0;381.0,12.0;382.0,12.0;383.0,12.0;384.0,12.0;385.0,12.0;386.0,12.0;387.0,12.0;388.0,12.0;389.0,12.0;390.0,12.0;391.0,12.0;392.0,12.0;393.0,12.0;394.0,12.0;395.0,12.0;396.0,12.0;397.0,12.0;398.0,12.0;399.0,12.0;400.0,12.0;401.0,12.0;402.0,12.0;403.0,12.0;404.0,12.0;405.0,12.0;406.0,12.0;407.0,12.0;408.0,12.0;409.0,12.0;410.0,12.0;411.0,12.0;412.0,12.0;413.0,12.0;414.0,12.0;415.0,12.0;416.0,12.0;417.0,12.0;418.0,12.0;419.0,12.0;420.0,12.0;421.0,12.0;422.0,12.0;423.0,12.0;424.0,12.0;425.0,12.0;426.0,12.0;427.0,12.0;428.0,12.0;429.0,12.0;430.0,12.0;431.0,12.0;432.0,12.0;433.0,12.0;434.0,12.0;435.0,12.0;436.0,12.0;437.0,12.0;438.0,12.0;439.0,12.0;440.0,12.0;441.0,12.0;442.0,12.0;443.0,12.0;444.0,12.0;445.0,12.0;446.0,12.0;447.0,12.0;448.0,12.0;449.0,12.0;450.0,12.0;451.0,12.0;452.0,12.0;453.0,12.0;454.0,12.0;455.0,12.0;456.0,12.0;457.0,12.0;458.0,12.0;459.0,12.0;460.0,12.0;461.0,12.0;462.0,12.0;463.0,12.0;464.0,12.0;465.0,12.0;466.0,12.0;467.0,12.0;468.0,12.0;469.0,12.0;470.0,12.0;471.0,12.0;472.0,12.0;473.0,12.0;474.0,12.0;475.0,12.0;476.0,12.0;477.0,12.0;478.0,12.0;479.0,12.0;480.0,12.0;481.0,12.0;482.0,12.0;483.0,12.0;484.0,12.0;485.0,12.0;486.0,12.0;487.0,12.0;488.0,12.0;489.0,12.0;490.0,12.0;491.0,12.0;492.0,12.0;493.0,12.0;494.0,12.0;495.0,12.0;496.0,12.0;497.0,12.0;498.0,12.0;499.0,12.0;500.0,12.0;501.0,12.0;502.0,12.0;503.0,12.0;504.0,12.0;505.0,12.0;506.0,12.0;507.0,12.0;508.0,12.0;509.0,12.0;510.0,12.0;511.0,12.0;512.0,12.0;513.0,12.0;514.0,12.0;515.0,12.0;516.0,12.0;517.0,12.0;518.0,12.0;519.0,12.0;520.0,12.0;521.0,12.0;522.0,12.0;523.0,12.0;524.0,12.0;525.0,12.0;526.0,12.0;527.0,12.0;528.0,12.0;529.0,12.0;530.0,12.0;531.0,12.0;532.0,12.0;533.0,12.0;534.0,12.0;535.0,12.0;536.0,12.0;537.0,12.0;538.0,12.0;539.0,12.0;540.0,12.0;541.0,12.0;542.0,12.0;543.0,12.0;544.0,12.0;545.0,12.0;546.0,12.0;547.0,12.0;548.0,12.0;549.0,12.0;550.0,12.0;551.0,12.0;552.0,12.0;553.0,12.0;554.0,12.0;555.0,12.0;556.0,12.0;557.0,12.0;558.0,12.0;559.0,12.0;560.0,12.0;561.0,12.0;562.0,12.0;563.0,12.0;564.0,12.0;565.0,12.0;566.0,12.0;567.0,12.0;568.0,12.0;569.0,12.0;570.0,12.0;571.0,12.0;572.0,12.0;573.0,12.0;574.0,12.0;575.0,12.0;576.0,12.0;577.0,12.0;578.0,12.0;579.0,12.0;580.0,12.0;581.0,12.0;582.0,12.0;583.0,12.0;584.0,12.0;585.0,12.0;586.0,12.0;587.0,12.0;588.0,12.0;589.0,12.0;590.0,12.0;591.0,12.0;592.0,12.0;593.0,12.0;594.0,12.0;595.0,12.0;596.0,12.0;597.0,12.0;598.0,12.0;599.0,12.0;600.0,12.0;601.0,12.0;602.0,12.0;603.0,12.0;604.0,12.0;605.0,12.0;606.0,12.0;607.0,12.0;608.0,12.0;609.0,12.0;610.0,12.0;611.0,12.0;612.0,12.0;613.0,12.0;614.0,12.0;615.0,12.0;616.0,12.0;617.0,12.0;618.0,12.0;619.0,12.0;620.0,12.0;621.0,12.0;622.0,12.0;623.0,12.0;624.0,12.0;625.0,12.0;626.0,12.0;627.0,12.0;628.0,12.0;629.0,12.0;630.0,12.0;631.0,12.0;632.0,12.0;633.0,12.0;634.0,12.0;635.0,12.0;636.0,12.0;637.0,12.0;638.0,12.0;639.0,12.0;640.0,12.0;641.0,12.0;642.0,12.0;643.0,12.0;644.0,12.0;645.0,12.0;646.0,12.0;647.0,12.0;648.0,12.0;649.0,12.0;650.0,12.0;651.0,12.0;652.0,12.0;653.0,12.0;654.0,12.0;655.0,12.0;656.0,12.0;657.0,12.0;658.0,12.0;659.0,12.0;660.0,12.0;661.0,12.0;662.0,12.0;663.0,12.0;664.0,12.0;665.0,12.0;666.0,12.0;667.0,12.0;668.0,12.0;669.0,12.0;670.0,12.0;671.0,12.0;672.0,12.0;673.0,12.0;674.0,12.0;675.0,12.0;676.0,12.0;677.0,12.0;678.0,12.0;679.0,12.0;680.0,12.0;681.0,12.0;682.0,12.0;683.0,12.0;684.0,12.0;685.0,12.0;686.0,12.0;687.0,12.0;688.0,12.0;689.0,12.0;690.0,12.0;691.0,12.0;692.0,12.0;693.0,12.0;694.0,12.0;695.0,12.0;696.0,12.0;697.0,12.0;698.0,12.0;699.0,12.0;700.0,12.0;701.0,12.0;702.0,12.0;703.0,12.0;704.0,12.0;705.0,12.0;706.0,12.0;707.0,12.0;708.0,12.0;709.0,12.0;710.0,12.0;711.0,12.0;712.0,12.0;713.0,12.0;714.0,12.0;715.0,12.0;716.0,12.0;717.0,12.0;718.0,12.0;719.0,12.0;720.0,12.0;721.0,12.0;722.0,12.0;723.0,12.0;724.0,12.0;725.0,12.0;726.0,12.0;727.0,12.0;728.0,12.0;729.0,12.0;730.0,12.0;731.0,12.0;732.0,12.0;733.0,12.0;734.0,12.0;735.0,12.0;736.0,12.0;737.0,12.0;738.0,12.0;739.0,12.0;740.0,12.0;741.0,12.0;742.0,12.0;743.0,12.0;744.0,12.0;745.0,12.0;746.0,12.0;747.0,12.0;748.0,12.0;749.0,12.0;750.0,12.0;751.0,12.0;752.0,12.0;753.0,12.0;754.0,12.0;755.0,12.0;756.0,12.0;757.0,12.0;758.0,12.0;759.0,12.0;760.0,12.0;761.0,12.0;762.0,12.0;763.0,12.0;764.0,12.0;765.0,12.0;766.0,12.0;767.0,12.0;768.0,12.0;769.0,12.0;770.0,12.0;771.0,12.0;772.0,12.0;773.0,12.0;774.0,12.0;775.0,12.0;776.0,12.0;777.0,12.0;778.0,12.0;779.0,12.0;780.0,12.0;781.0,12.0;782.0,12.0;783.0,12.0;784.0,12.0;785.0,12.0;786.0,12.0;787.0,12.0;788.0,12.0;789.0,12.0;790.0,12.0;791.0,12.0;792.0,12.0;793.0,12.0;794.0,12.0;795.0,12.0;796.0,12.0;797.0,12.0;798.0,12.0;799.0,12.0;800.0,12.0;801.0,12.0;802.0,12.0;803.0,12.0;804.0,12.0;805.0,12.0;806.0,12.0;807.0,12.0;808.0,12.0;809.0,12.0;810.0,12.0;811.0,12.0;812.0,12.0;813.0,12.0;814.0,12.0;815.0,12.0;816.0,12.0;817.0,12.0;818.0,12.0;819.0,12.0;820.0,12.0;821.0,12.0;822.0,12.0;823.0,12.0;824.0,12.0;825.0,12.0;826.0,12.0;827.0,12.0;828.0,12.0;829.0,12.0;830.0,12.0;831.0,12.0;832.0,12.0;833.0,12.0;834.0,12.0;835.0,12.0;836.0,12.0;837.0,12.0;838.0,12.0;839.0,12.0;840.0,12.0;841.0,12.0;842.0,12.0;843.0,12.0;844.0,12.0;845.0,12.0;846.0,12.0;847.0,12.0;848.0,12.0;849.0,12.0;850.0,12.0;851.0,12.0;852.0,12.0;853.0,12.0;854.0,12.0;855.0,12.0;856.0,12.0;857.0,12.0;858.0,12.0;859.0,12.0;860.0,12.0;861.0,12.0;862.0,12.0;863.0,12.0;864.0,12.0;865.0,12.0;866.0,12.0;867.0,12.0;868.0,12.0;869.0,12.0;870.0,12.0;871.0,12.0;872.0,12.0;873.0,12.0;874.0,12.0;875.0,12.0;876.0,12.0;877.0,12.0;878.0,12.0;879.0,12.0;880.0,12.0;881.0,12.0;882.0,12.0;883.0,12.0;884.0,12.0;885.0,12.0;886.0,12.0;887.0,12.0;888.0,12.0;889.0,12.0;890.0,12.0;891.0,12.0;892.0,12.0;893.0,12.0;894.0,12.0;895.0,12.0;896.0,12.0;897.0,12.0;898.0,12.0;899.0,12.0;900.0,12.0;901.0,12.0;902.0,12.0;903.0,12.0;904.0,12.0;905.0,12.0;906.0,12.0;907.0,12.0;908.0,12.0;909.0,12.0;910.0,12.0;911.0,12.0;912.0,12.0;913.0,12.0;914.0,12.0;915.0,12.0;916.0,12.0;917.0,12.0;918.0,12.0;919.0,12.0;920.0,12.0;921.0,12.0;922.0,12.0;923.0,12.0;924.0,12.0;925.0,12.0;926.0,12.0;927.0,12.0;928.0,12.0;929.0,12.0;930.0,12.0

Q. Which is the most informative pattern?
Q. Which technique provides it?




Fingerprinting on most relevant chemical classes - Aldehydes



Homologue series fingerprinting – structured patterns

TIC
D-HS-PDMS
Headspace




Odor thresholds


Undecanal: 5 µg/L

2-Undecanone: 7 µg/L

(E)-2-Decenal: 0.3-0.4 µg/L




TIC
MMSE
Adsorption
In-solution





Fingerprinting on most relevant chemical classes


Lactones




TIC
HSSE-PDMS



TIC zoom
HSSE-PDMS



SIM
HSSE-PDMS



Q. Which are analytes that represent the Chemical Odor Code of milk?

Modern omics disciplines
Advanced Sensomics




4D Odor

1D-2D-3D Chemical Code

ORs Code

Odor descriptor

- cabbage
- cheesy
- fatty
- green
- hay-cucumber
- rancid
- rubber
- sour-vinegar
- stale, hay
- sweaty
- sweet
- sweet-fatty
- sweet-green
- vitamin-rubber
- waxy
- waxy-green




- decanoic acid
- 2-furanmethanol
- acetic acid
- maltol
- octanoic acid
- decenoic acid
- (E,E)-2,4-decadienal
- butanoic acid
- dodecanoic acid
- hexanal
- (E)-2-undecenal
- (E)-2-nonenal
- 3-methylbutanal
- hexanoic acid
- 3-methylbutanoic acid
- nonanoic acid
- β-ionone
- pentanoic acid
- benzothiazole
- dimethyltrisulfide
- δ-undecalactone
- δ-decalactone
- γ-dodecalactone


- CeCR472
- ORL1138
- ORL1500
- ORL1520
- ORL1549
- ORL1580
- ORL1582
- ORL1629
- ORL2156
- ORL2162
- ORL3370
- ORL461
- ORL464
- ORL465
- ORL467
- ORL468
- ORL470
- ORL471
- ORL472
- ORL52
- ORL533
- ORL876

Experimental data

Prediction




 **Advanced *Sensomics*** and its peculiar investigation strategies are **fundamental to identify triggering factors** (ligands-odorants) **of odor perception**


 Our contribution rely on the **ability to consider all possible informative dimensions** and **giving** them the **appropriate priority** during investigation

 Advanced analytical **complex** platforms could effectively contribute to unlock closed doors giving **new perspectives in biological phenomena interpretation**

Odorants play also a role in human metabolism regulation
Behind pleasure we can find a further key for well-being



Results of bio-activation of three flavor active compounds on the PPAR γ receptor.




Audouze K, Tromelin A, Le Bon AM, Belloir C, et al. (2014) Identification of Odorant-Receptor Interactions by Global Mapping of the Human Odorome. PLoS ONE 9(4): e93037. doi:10.1371/journal.pone.0093037

<http://www.plosone.org/article/info:doi/10.1371/journal.pone.0093037>



Grazie per l'attenzione



*Carlo
Bicchi*



*Patrizia
Rubiolo*



*Cecilia
Cagliero*



*Barbara
Sgorbini*



*Erica
Liberto*

